

***InVitroGRO*<sup>™</sup> PT Medium**

Product No.	Description	Size
Z990001	<i>InVitroGRO</i> PT Medium	250 mL

**Product Description**

*InVitroGRO* PT Medium is formulated to facilitate the optimal growth of In Vitro Technologies human renal proximal tubule cells. *InVitroGRO* PT Medium is a complete medium that does not contain antibiotics. It is recommended that antibiotics be added prior to use.

**Stability:** 60 days for *InVitroGRO* PT Medium that does not contain antibiotics.

**Storage:** 4°C

**Caution:** Treat all products containing human and monkey-derived materials a potentially infectious, as no known test methods can offer assurance that products derived from human or monkey tissues will not transmit infectious agents.

All products are for research use only. Do not use in animals or humans. These products have not been approved for any diagnostic or clinical procedures.