

***InVitroGRO*[™] HI Medium**

Product No.	Description	Size
Z99009	<i>InVitroGRO</i> HI Medium	250 mL
Z990012	<i>InVitroGRO</i> HI Medium	500 mL
Z990002	<i>InVitroGRO</i> HI Medium	1 L

Product Description

InVitroGRO HI Medium is a hepatocyte incubation medium used for culturing primary hepatocytes¹. *InVitroGRO* HI Medium is formulated without antibiotics to increase its shelf life. In Vitro Technologies recommends the use of the *Torpedo*[™] Antibiotic Mix with *InVitroGRO* HI Medium.

Stability: 90 days for *InVitroGRO* HI Medium without *Torpedo* Antibiotic Mix
7 days for Complete *InVitroGRO* HI Medium

Storage: 4°C

Reference

1. Roymans, D.; Van Looveren, C.; Leone, A.; Parker, J. B.; McMillan, M.; Johnson, M. D.; Koganti, A.; Gilissen, R.; Silber, P.; Mannens, G.; Meuldermans, W. Determination of cytochrome P450 1A2 and cytochrome P450 3A4 induction in cryopreserved human hepatocytes. *Biochem. Pharmacol.* **2004**, *67*(3), 427–437.

Caution: Treat all products containing human and monkey-derived materials a potentially infectious, as no known test methods can offer assurance that products derived from human or monkey tissues will not transmit infectious agents.

All products are for research use only. Do not use in animals or humans. These products have not been approved for any diagnostic or clinical procedures.