

**Fact Sheet 303*****“Celsis RapiScreen™ Biologics Kit Detects a Wide Range of Organisms\*”***

Celsis RapiScreen™ Biologics kit can detect a wide range of organisms relevant to the bioprocessing industry. Micro-organisms examined include:

*Bacillus subtilis* ATCC 6633

*Burkholderia cepacia* ATCC 25416

*Candida albicans* ATCC 10231

*Enterococcus faecium* ATCC 35667

*Enterococcus gallinarum* ATCC 700425

*Escherichia coli* ATCC 9739

*Pseudomonas aeruginosa* ATCC 9027

*Ralstonia pickettii* ATCC 49129

*Staphylococcus aureus* ATCC 6538

*Staphylococcus epidermidis* ATCC 12228

\*Organisms detected in cell suspensions added to enrichment broth.